



Contents

1	Introduction	2
2	Prerequisites	2
2.1	Required Python Modules	3
2.2	Register ZEN Functionality	4
3	ZEN Blue - Setting up the Experiment	5
3.1	ZEN - Configure a Sample Camera	6
4	Python - Control ZEN from a python script	7
4.1	Control ZEN from a Python script	8
4.2	ZEN-Python at work	9
4.3	Importing CZI files into Python	11
4.4	Importing the CZI image data into a NumPy array	11
4.5	Displaying the analysis results	11
5	Appendix: Python Scripts	12
5.1	RunZenfromPython.py	12
5.2	Additional Python Scripts	13
5.2.1	Package BioFormatsRead	13
5.2.2	ReadAnalyzeImageData.py	14
5.2.3	wellplatetools.py	16
6	Acknowledgements and Links	17
7	Disclaimer	18



1 Introduction

This application note will explain how to create a workflow using the ZEN-Python connection. The basic idea is to control everything from within Python (Master). While ZEN internally is using IronPython, this notes describes how to control ZEN from within a CPython implementation

The ZEN Image Acquisition (Slave) software is "only" doing the image acquisition. The signal to start the experiment is send from Python to ZEN.

When the experiment is finished, the CZI data are imported into Python using BioFormats and "some" simple image analysis is carried out to underline the workflow concept.

2 Prerequisites

To be able to use the functionality described herein, you need

- **ZEN Blue 2012 or later with the Macro Environment module** - http://www.zeiss.com/microscopy/en_us/products/microscope-software/zen-2012.html
- **Python 3.5 64-bit Installation with some extensions** - <https://www.python.org/>

To profit fully from the following example, you should have working experience with both applications. You should be familiar with the user interface of your favorite python IDE and the creation and use of scripts.



2.1 Required Python Modules

In order to try out the complete example including the image analysis part some additional python modules are required. Most of the packages are included into typical python distribution like Anaconda (<https://www.continuum.io/downloads>). Some of the additional other modules, that might not be part of the python distribution can be typically found on PyPI (<https://pypi.python.org/pypi>).

- **python-bioformats**
- **Javabridge**
- **tiffle**
- **mahotas**
- **czifile.py** - <http://www.lfd.uci.edu/~gohlke>



2.2 Register ZEN Functionality

To be able to use ZEN services in a COM environment, the ZEN functionality must be registered as follows using the following BAT-File, which **must be executed as administrator**.

```
1  :: regScripting_Release.bat
2
3  echo off
4
5  pushd "C:\Windows\Microsoft.NET\Framework64\v4.0.30319"
6
7  SET dll-1="C:\Program Files\Carl Zeiss\ZEN 2\ZEN 2 (blue edition)\Zeiss.Micro.Scripting.dll"
8  regasm /u /codebase /tlb %dll-1%
9  regasm /codebase /tlb %dll-1%
10
11 SET dll-2="C:\Program Files\Carl Zeiss\ZEN 2\ZEN 2 (blue edition)\Zeiss.Micro.LM.Scripting.dll"
12 regasm /u /codebase /tlb %dll-2%
13 regasm /codebase /tlb %dll-2%
14 popd
15 pause
```

The final result of this set of commands is that the .NET classes and structures within **Zeiss.Micro.Scripting.dll** and **Zeiss.Micro.LM.Scripting.dll** are made known and accessible within the COM environment as well.

Remark: Please note that you must edit the BAT-file and adapt the path location, where the DLLs can be found to your needs.

3 ZEN Blue - Setting up the Experiment

The easiest way to setup a wellplate experiment is to use the Tile & Position Module. It is assumed that one is familiar with setting up ZEN Blue experiments in general since this is not explained in detail here.

- Specify the correct folder for your sample camera in case of a simulated experiment.
- Use Smart Setup to create a (DAPI) channel.
- Activate **Tiles**, change to **Position Setup** and select **Carrier** mode.
- Add **one** position per well.

The ZEN Blue software should look similar to this:




Figure 1: ZEN Blue - Wellplate Experiment Setup

3.1 ZEN - Configure a Sample Camera

The idea is to place some "real sample" images in a specific folder in order to use them for a simulated experiment including a real image analysis. This is an extremely useful feature to test and demonstrate most of the capabilities of ZEN Blue Experiments without the need to have real hardware in place.

This requires a special XML-file which is usually located at C:/ProgramData/Carl Zeiss/MTB2011/2.7.0.5/CZIS_Cameras.xml (Windows 7, depending on the current ZEN Blue and MTB version).

If configured correctly, one can specify the image folder inside ZEN on the Camera Control on the Locate tab or the Acquisition Mode ToolWindow on the Acquisition Tab:




Figure 2: ZEN Blue - Sample Camera Configuration



4 Python - Control ZEN from a python script

The main idea behind controlling ZEN from within Python is to make Python aware of the ZEN Scripting DLLs via the .COM interface. This way you can use all the commands available inside the OAD python scripts inside ZEN directly within a normal python script inside your favorite IDE.

The main ToDo's in Python inside this application note are:

- **Import the ZEN objects into Python and run the actual experiment.**
- **Read the CZI file using BioFormats into a numpy array.**
- **Cycle through all the images and do some "meaningful" analysis.**
- **Display the image analysis results.**

Of course the most important part of this example is to understand the implications of the 1st point. Once the ZEN-Python bridge is established, the ZEN OAD Simple-API Interface becomes part of the running python script and can be used in conjunction with the rest of python functionality.

4.1 Control ZEN from a Python script

This is the main python script controlling the complete workflow. It establishes the connection to ZEN by making the OAD Simple-API available within python.

```
1  # -*- coding: utf-8 -*-
2  """
3  File: RunZenfromPython.py
4  Date: 03.05.2017
5  Author: Sebastian Rhode
6  Version: 1.1
7
8  This Python script demonstrates the capabilities of the .COM interface
9  to establish a connection between ZEN Blue and Python.
10 This connection allows to use ZEN Blue OAD Simple-API objects from
11 within a python script.
12
13 Requirements:
14 - bftools.py, czitools.py, miscctool.py
15 - tiffiofile
16 - czifile.py
17 - mahotas
18 """
19
20
21 import win32com.client
22 import os
23 import ReadAnalyzeImageData as rad
24
25 # Define place to store the CZI file
26 savefolder = 'C:\\Python_ZEN_Output\\'
27 # check if the folder already exists
28 try:
29     os.makedirs(savefolder)
30 except OSError:
31     if not os.path.isdir(savefolder):
32         raise
33
34 # Import the ZEN OAD Scripting into Python
35 Zen = win32com.client.GetActiveObject("Zeiss.Micro.Scripting.ZenWrapperLM")
```

The most important line is number 35. Here you tell Python to load the ZEN Scripting API. Once this is done, you can use all available ZEN objects (from the Simple API) the same way one would do it inside the ZEN software.

- Define save folder and ZEN Experiment.
- Get the ZEN experiment using the correct name and run it.
- Read the CZI file using python-bioformats.
- Analyze the image from every well and display the results.

Of course there is little bit more behind the scenes than shown inside this script, but all the functionality to control ZEN from within Python is located inside this script.

This opens up a whole range of exciting uses cases. You can use all the functionality of Python for image and data analysis or even hardware control while having the option to use the powerful ZEN Blue acquisition engine for all your imaging experiments and controlling the microscope itself.

4.2 ZEN-Python at work

Now we are ready to start the workflow by running **RunZenfromPython.py** from within the python IDE and watch the output.




Figure 3: ZEN and Python running

Now one just has to start the python script inside the IDE. This will launch the actual ZEN Experiment and finally read and analyze the resulting CZI file in order to display some plots.

```

IPython console
Console 2/A
Python 3.5.2 [Anaconda custom (64-bit)] (default, Jul 5 2016, 11:41:13) [MSC v.1900 64 bit (AMD64)]
Type "copyright", "credits" or "license()" for more information.

IPython 5.1.0 -- An enhanced Interactive Python.
? -> Introduction and overview of IPython's features.
%quickref -> Quick reference.
help -> Python's own help system.
object? -> Details about 'object', use 'object??' for extra details.

In [1]: runfile('C:/TFS/Doc/3-ZIS/2-ProjectManagement/Teilprojekte/Teilprojekt_Dokumentation/Temporary Doku (raw
data)/Rhode/Python/Appnote_Control_ZEN_Python/SourceCode/RunZenfromPython.py', wdir='C:/TFS/Doc/3-ZIS/2-
ProjectManagement/Teilprojekte/Teilprojekt_Dokumentation/Temporary Doku (raw data)/Rhode/Python/Appnote_Control_ZEN_Python/SourceCode')
Multi-Resolution API not enabled yet.
Retrieving Image Dimensions ...
T: 1 Z: 1 C: 1 X: 640 Y: 640
Using czifile.py to get CZI Shape info.
Using BioFormats to get MetaInformation.

OME NameSpace used : http://www.openmicroscopy.org/Schemas/OME/2016-06
BF Version used : c:\Users\MISRH\Documents\Software\BioFormats_Package\5.4.1\bioformats_package.jar
-----
Image Directory : C:\Python_ZEN_Output
Image Filename : Experiment-36.czi
MultiResolution : False
Images Dim Sizes [0] : [96, 1, 1, 1, 640, 640]
Dimension Order BF : XYZCT
Dimension Order CZI : b'85CVX0'
Shape CZI : (1, 96, 1, 158140, 248140, 1)
Total Series Number : 96
Image Dimensions : 96 1 1 1 640 640
Scaling XYZ [micron] : 0.4 0.4 0.4
Objective M-NA-Imm : 5.0 0.35 Air
Objective Name : Plan-Apochromat 5x/0.35
Ex. Wavelengths [nm] : [353.]
Em. Wavelengths [nm] : [485.]
Dyes : ['None']
Detector Model : TestCam
Detector Name : Detector:Internal
Detector ID : Detector:Internal
Channels : ['DAPI']
Channel Description : []
ImageIDs : [0, 1, 2, 3, 4, 5, 6, 7, 8, 9, 10, 11, 12, 13, 14, 15, 16, 17, 18, 19, 20, 21, 22, 23, 24, 25, 26, 27, 28, 29, 30, 31,
32, 33, 34, 35, 36, 37, 38, 39, 40, 41, 42, 43, 44, 45, 46, 47, 48, 49, 50, 51, 52, 53, 54, 55, 56, 57, 58, 59, 60, 61, 62, 63, 64, 65, 66,
67, 68, 69, 70, 71, 72, 73, 74, 75, 76, 77, 78, 79, 80, 81, 82, 83, 84, 85, 86, 87, 88, 89, 90, 91, 92, 93, 94, 95]
Start Processing ...
Done!

```

Figure 4: Python - Output of the running script inside console




Figure 5: Python - Visualization of Image Analysis Results



4.3 Importing CZI files into Python

ZEN Blue uses the CZI data format to store the data. Due to its open concept it is very well supported by the BioFormats library. Further info can be found at:

- www.zeiss.com/czi

Instructions on the usage of the Python wrapper for BioFormats can be found here:

- <http://pythonhosted.org/python-bioformats>
- <https://github.com/CellProfiler/python-bioformats/>

There is more than one way to import a CZI file, so the choice depends on the nature of the data. For this example we use a pretty generic approach to read the CZI dataset.

4.4 Importing the CZI image data into a NumPy array

A very important topic are the CZI Metadata. Again we can rely on the BioFormats library and just use the already existing functionality to get all the information we need. For the example we use the `get_image6d` function provided by the `bftools.py` script.

4.5 Displaying the analysis results

To make this example more descriptive, a few additional scripts will be used to display the data and demonstrate the idea to automate a complete workflow within Python. ZEN is a very powerful image acquisition software, but especially when it comes data analysis and crunching numbers, one may needs to leave the ZEN world. For this example the python package `mahotas` was used to do the image analysis and `matplotlib` was used to create the plots. The python script to create the simple heatmap can be found in **5.2.3**.

5 Appendix: Python Scripts

5.1 RunZenfromPython.py

This is the main Python script controlling the complete workflow.

```
1  # -*- coding: utf-8 -*-
2  """
3  File: RunZenfromPython.py
4  Date: 03.05.2017
5  Author: Sebastian Rhode
6  Version: 1.1
7
8  This Python script demonstrates the capabilities of the .COM interface
9  to establish a connection between ZEN Blue and Python.
10 This connection allows to use ZEN Blue OAD Simple-API objects from
11 within a python script.
12
13 Requirements:
14
15 - bftools.py, czitools.py, miscctool.py
16 - tiffiofile
17 - czifile.py
18 - mahotas
19 """
20
21 import win32com.client
22 import os
23 import ReadAnalyzeImageData as rad
24
25 # Define place to store the CZI file
26 savefolder = 'C:\\Python_ZEN_Output\\'
27 # check if the folder already exists
28 try:
29     os.makedirs(savefolder)
30 except OSError:
31     if not os.path.isdir(savefolder):
32         raise
33
34 # Import the ZEN OAD Scripting into Python
35 Zen = win32com.client.GetActiveObject("Zeiss.Micro.Scripting.ZenWrapperLM")
36
37 # Define the experiment to be executed
38 ZEN_Experiment = "ML_96_Wellplate_Castor.czexp"
39
40 # run the experiment in ZEN and save the data to the specified folder
41 exp = Zen.Acquisition.Experiments.GetByName(ZEN_Experiment)
42 img = Zen.Acquisition.Execute(exp)
43
44 # Show the image in ZEN
45 Zen.Application.Documents.Add(img)
46
47 # Use the correct save method - it is polymorphic ... :)
48 filename = savefolder + img.Name_2
49 img.Save_2(filename)
50
51 # get the actual CZI image data using Python wrapper for BioFormats
52 img6d = rad.ReadImage(filename)
53
54 # Analyze the images - Example: Count Cells
55 obj, labeled = rad.CountObjects(img6d)
56
57 # Display some data
58 rad.DisplayData_2(obj, labeled)
```



5.2 Additional Python Scripts

This main script uses the following additional tools:

- `bftools.py`
- `czitools.py`
- `bftools.py`
- `ReadAnalyseImageData.py`
- `wellplatetool.py`

5.2.1 Package BioFormatsRead

This is a little collection of python tools used to actually read the CZI image data via BioFormats. It used the python wrapper of the BioFormats library and `czifile.py` to read CZI image data into a NumPy array.

The scripts and additional tools can be also found here:

<https://github.com/sebi06/BioFormatsRead>

For more detailed information on the usage please check the readme on the project page

5.2.2 ReadAnalyzeImageData.py

```

1  # -*- coding: utf-8 -*-
2  """
3  File: ReadAnalyzeImageData.py
4  Date: 03.05.2017
5  Author: Sebastian Rhode
6  Version: 0.3
7  """
8
9  from matplotlib import pyplot as plt, cm
10 import numpy as np
11 import mahotas
12 import wellplatetool as wp
13 import sys
14 import bftools as bf
15
16
17 def ReadImage(filename):
18
19     urlnamespace = 'http://www.openmicroscopy.org/Schemas/OME/2016-06'
20     # specify bioformats_package.jar to use if required
21     bfpackage = r'c:\Users\MISRH\Documents\Software\BioFormats_Package\5.4.1\bioformats_package.jar'
22     bf.set_bfpath(bfpackage)
23
24     # get image meta-information
25     MetaInfo = bf.get_relevant_metainfo_wrapper(filename, namespace=urlnamespace, bfpath=bfpackage, showinfo=False)
26     img6d, readstate = bf.get_image6d(filename, MetaInfo['Sizes'])
27
28     # show relevant image Meta-Information
29     bf.showtypicalmetadata(MetaInfo, namespace=urlnamespace, bfpath=bfpackage)
30
31     return img6d
32
33
34 def CountObjects(img6d):
35
36     obj = np.zeros(img6d.shape[0])
37
38     print('Start Processing ...')
39     steps = img6d.shape[0]/10
40
41     # count cells with individual thresholds per frame
42     for i in range(0, img6d.shape[0], 1):
43
44         img = img6d[i, 0, 0, 0, :, :]
45         T = mahotas.otsu(img)
46         img = (img > T)
47         img = mahotas.gaussian_filter(img, 0.5)
48         labeled, numobjects = mahotas.label(img)
49         obj[i] = numobjects
50         if i%steps == 0:
51             print('\b.',)
52             sys.stdout.flush()
53
54     print('Done!')
55
56     return obj, labeled
57
58
59 def DisplayData_1(obj, labeled):
60
61     fig1 = plt.figure()
62     ax1 = fig1.add_subplot(111)
63     cax = ax1.matshow(labeled, interpolation='nearest', cmap=cm.jet)
64     cbar = fig1.colorbar(cax)
65     ax1.set_title('Example - Labeled Objects')
66     ax1.set_xlabel('X-dimension [pixel]')
67     ax1.set_ylabel('Y-dimension [pixel]')
68
69     fig2 = plt.figure()
70     ax2 = fig2.add_subplot(111)
71     ax2.hist(obj, normed=0, facecolor='green', alpha=0.75)
72     ax2.set_title('Distribution - Number of Objects')
73     ax2.set_xlabel('Number')
74     ax2.set_ylabel('Occurrence')
75     ax2.grid(True)
76
77     wp.ShowPlateData(obj, 8, 12, 'Object Number')
78
79     # show plots
80     plt.show()
81
82
83 def DisplayData_2(obj, labeled):
84

```

```

85 fig = plt.figure(figsize=(12, 8), dpi=100)
86 ax1 = fig.add_subplot(221)
87 ax1.hist(obj, normed=0, facecolor='green', alpha=0.75)
88 ax1.set_title('Distribution - Number of Objects')
89 ax1.set_xlabel('Number')
90 ax1.set_ylabel('Occurrence')
91 ax1.grid(True)
92
93 ax2 = fig.add_subplot(222)
94 ind = np.arange(obj.shape[0])
95 ax2.bar(ind, obj, width=0.8, color='r')
96 ax2.set_title('Follow Acquisition')
97 ax2.set_ylabel('Number')
98 ax2.set_xlim(0, obj.shape[0])
99 ax2.grid(True)
100
101 ax3 = fig.add_subplot(223)
102 Nc = 12
103 Nr = 8
104 WellData, labelx, labely, fs = wp.ReturnPlateHeatmap(obj, Nr, Nc)
105 cax = ax3.imshow(WellData, interpolation='nearest', cmap=cm.jet)
106 cbar = fig.colorbar(cax)
107 ax3.set_title('Objects per Well')
108 ax3.set_xticks(np.arange(0, Nc, 1))
109 ax3.set_xticklabels(labelx, fontsize=fs)
110 ax3.set_yticks(np.arange(0, Nr, 1))
111 ax3.set_yticklabels(labely, fontsize=fs)
112
113 ax4 = fig.add_subplot(224)
114 cax = ax4.matshow(labeled, interpolation='nearest', cmap=cm.jet)
115 ax4.set_xlabel('X-dimension [pixel]')
116 ax4.set_ylabel('Y-dimension [pixel]')
117
118 # show plots
119 plt.show()

```

5.2.3 wellplatetools.py

```

1  # -*- coding: utf-8 -*-
2  """
3  File: wellplatetool.py
4  Date: 03.05.2017
5  Author: Sebastian Rhode
6  Version: 1.2
7  """
8
9  import numpy as np
10 import matplotlib.pyplot as plt
11 from matplotlib import cm
12
13
14 def ExtractLabels(Nr, Nc):
15
16     # labeling schemes
17     LabelX = ['01','02','03','04','05','06','07','08','09','10','11','12',
18              '13','14','15','16','17','18','19','20','21','22','23','24',
19              '25','26','27','28','29','30','31','32','33','34','35','36',
20              '37','38','39','40','41','42','43','44','45','46','47','48',]
21
22     LabelY = ['A','B','C','D','E','F','G','H','I','J','K','L','M','N','O','P',
23              'Q','R','S','T','U','V','W','X','Y','Z','AA','AB','AC','AD','AE','AF']
24
25     Lx = LabelX[0:Nc]
26     Ly = LabelY[0:Nr]
27
28     return Lx, Ly
29
30
31 def ShowPlateData(data, Nr, Nc, parameter):
32
33     # reshape data to create 2D matrix
34     WellData = data.reshape(Nr, Nc)
35     [labelx, labely] = ExtractLabels(Nr, Nc)
36
37     # create figure
38     fig = plt.figure(figsize=(6, 4), dpi=100)
39     ax1 = fig.add_subplot(111)
40     # set colormap
41     cmap = cm.jet
42     # show the well plate as an image
43     cax = ax1.imshow(WellData, interpolation='nearest', cmap=cmap)
44     # determine an appropriate font size
45     if Nr <= 32 and Nr > 16:
46         fs = 7
47     elif Nr <= 16 and Nr > 8:
48         fs = 9
49     elif Nr <= 8:
50         fs = 11
51     # format the display
52     ax1.set_xticks(np.arange(0, Nc, 1))
53     ax1.set_xticklabels(labelx, fontsize=fs)
54     ax1.set_yticks(np.arange(0, Nr, 1))
55     ax1.set_yticklabels(labely, fontsize=fs)
56     ax1.set_title(parameter)
57     cbar = fig.colorbar(cax)
58
59
60 def ReturnPlateHeatmap(data, Nr, Nc):
61
62     # reshape data to create 2D matrix
63     WellData = data.reshape(Nr, Nc)
64     [labelx, labely] = ExtractLabels(Nr, Nc)
65
66     ## determine an appropriate font size
67     if Nr <= 32 and Nr > 16:
68         fontsize = 7
69     elif Nr <= 16 and Nr > 8:
70         fontsize = 9
71     elif Nr <= 8:
72         fontsize = 11
73
74     return WellData, labelx, labely, fontsize

```


6 Acknowledgements and Links

Special thanks to Christoph Gohlke for providing the Windows binaries for the providing the python extension packages at:

- <http://www.lfd.uci.edu/~gohlke/pythonlibs/>

and the CellProfiler and the BioFormats Team for providing the BioFormats library and the respective python wrapper.

- <https://github.com/CellProfiler/python-bioformats/>
- [www.http://loci.wisc.edu/software/bio-formats](http://loci.wisc.edu/software/bio-formats)

Thanks to Luis Pedro Coelho for the mahotas package.

Parts of the **bftools.py** script are based on the project lesion by Juan Nunez-Iglesias:

- <https://github.com/jni/lesion>



7 Disclaimer

This is an application note free to use for everybody. Use it on your own risk.

Carl Zeiss Microscopy GmbH's ZEN software allows to connect to a the third party software Python. Therefore Carl Zeiss Microscopy GmbH undertakes no warranty concerning Python, makes no representation that Python will work on your hardware and will not be liable for any damages caused by the use of this extension. By running this example you agree to this disclaimer.